

Imagen® Fast SYBR Green qPCR Premix

Cat. NO IPQ

Description


Imagen® Fast SYBR Green qPCR Premix is a ready-to-use, 2X concentrated premix reagent including a novel Hotstart Taq DNA polymerase, SYBR®Green I, optimized reaction buffer and nucleotides for running quantitative real-time PCR, including qPCR and 2-step qRT-PCR, in the SYBR®Green I detection format. Fast SYBR Green PCR Premix provides high sensitivity, wide dynamic range and reproducibility, quantification.

Contents

The **Imagen® Fast SYBR Green qPCR Premix** is supplied as a ready-to-use 2x reaction mix. The formulation contains Fast Hotstart Taq DNA polymerase, MgCl₂, dNTPs, reaction enhancers, and stabilizers.

Storage

- ✓ -20 °C
- ✓ Protected from light
- ✓ Avoid repeated freezing and thawing



Reaction Mix Thawing and Handling

Imagen® Fast SYBR Green qPCR Premix is delivered in a 2x ready-to-use format. To use the mix, thaw the vial on ice to 4 °C.

Please completely mix the vial and briefly centrifuge to ensure all components are at the bottom of the tube. Store on ice protected from light until ready to use. If using automated liquid handling, let sit at ambient temperature for 10 min to further reduce the viscosity.

Application

- ✓ Quantitative real-time PCR
- ✓ Quantitative 2-step RT-PCR
- ✓ Gene expression analysis
- ✓ Low copy gene detection

Prepare the qPCR Reaction Mix

1. Mix the **Imagen® Fast SYBR Green qPCR Premix** thoroughly but gently until it's completely homogenous.
2. Prepare the qPCR Reaction Mix for the number of reactions required as shown in table below and plus 10% overage.

Reagent	Volume (ul)	Final conc.
Imagen® Fast SYBR Green qPCR Premix	12.5	1x
Forward Primer(10 uM)	0.625	200 - 500 nM
Reverse Primer(10 uM)	0.625	200 - 500 nM
DNA Template	2	100 ng - 1 pg
Nuclease-free water	9.25	-
Final volume	25	-

3. Vortex the tube to mix the contents thoroughly, then centrifuge briefly to collect the contents at the bottom of the tube. (*Use good pipetting practice to ensure assay precision and accuracy of dispensing.)

4. Add DNA (and nuclease-free water, if needed) to the PCR tubes or wells containing the reaction mix, seal tubes or wells with flat caps or optically transparent film, and gently vortex to ensure thorough mixing of the reaction components.
5. Program the thermal cycling protocol on the real-time PCR instrument.

Step		Temp. °C	Time	Cycles
DNA polymerase activation and template denaturation		95°C	60 sec	1
Amplification	Template denaturation	95°C	3 sec	35-40
	Annealing / Extension and plate read	58 - 62°C	≥20 sec © Data acquisition	
Melt Curve		95°C 60°C 95°C	20 sec 30 sec 15 sec	1

6. Load the PCR tubes or plates onto the real-time PCR instrument and start the qPCR run program.
7. When thermal cycling is complete, perform data according to the instructions in the instrument-specific software.