

Imagen® TaqMan RT-qPCR Premix

Cat. No

Description


Imagen® TaqMan RT-qPCR Premix uses a designed fluorogenic probe to detect target PCR products during real-time PCR cycle, fluorescent is detectable when the quencher of probe is removed during the PCR extension cycle. Real-time TaqMan RT-PCR kit complete reverse transcription and TaqMan Real-time PCR in a single tube. Components necessary for performing TaqMan Real-time PCR from RNA template are specially prepared to ensure high specificity and high sensitivity.

Contents

The **Imagen® TaqMan RT-qPCR Premix** is supplied as a ready-to-use 2x reaction mix. The formulation contains, Hotstart Taq DNA polymerase, Reverse transcriptases, dNTPs, MgCl₂, reaction enhancers, and stabilizers.

Storage

- ✓ -20 °C
- ✓ Protected from light
- ✓ Avoid repeated freezing and throwing



Reaction Mix Thawing and Handling

Imagen® TaqMan RT-qPCR Premix is delivered in a 2x ready-to-use format. To use the mix, thaw the vial on ice to 4 °C.

Please completely mix the vial and briefly centrifuge to ensure all components are at the bottom of the tube. Store on ice protected from light until ready to use. If using automated liquid handling, let sit at ambient temperature for 10 min to further reduce the viscosity.

Application

- ✓ TaqMan RT-qPCR based on Specific Probes.
- ✓ Detection and Quantification of RNA targets.
- ✓ High Throughput Applications.

Prepare the qPCR Reaction Mix

1. Mix the **Imagen® TaqMan RT-qPCR Premix** thoroughly but gently until it's completely homogenous.
2. Prepare the qPCR Reaction Mix for the number of reactions required as shown in table below and plus 10% overage.
4. Add DNA (and nuclease-free water, if needed) to the PCR tubes or wells containing the reaction mix, seal tubes or wells with flat caps or optically transparent film, and gently vortex to ensure thorough mixing of the reaction components.

Reagent	Volume (ul)	Final conc.
Imagen® TaqMan RT-qPCR Premix	12.5	1x
Forward Primer(10 uM)	0.75	300 - 600 nM
Reverse Primer(10 uM)	0.75	300 - 600 nM
Fluorogenic Probe(10 uM)	0.5	200-400 nM
RNA Template	2	100 ng - 10 pg
Nuclease-free water	8.5	-
Final volume	25	-

3. Vortex the tube to mix the contents thoroughly, then centrifuge briefly to collect the contents at the bottom of the tube. (*Use good pipetting practice to ensure assay precision and accuracy of dispensing.)
5. Program the thermal cycling protocol on the real-time PCR instrument.
6. Load the PCR tubes or plates onto the real-time PCR instrument and start the qPCR run program.
7. When thermal cycling is complete, perform data according to the instructions in the instrument-specific software.

Step	Temp. °C	Time	Cycles
Reverse Transcription	48°C	30 min	1
DNA polymerase activation and template denaturation	95°C	10 min	1
Amplification	Template denaturation	95°C	20 sec
	Annealing / Extension and plate read	58 - 62°C	60 sec © Data acquisition